



**PERIPHERAL BLOOD FILM PREPARATION**

**1. Introduction**

Preparation of blood film is a fundamental technique in Hematology and blood collection is an important pre-analytical component. Red cells morphology and differential leukocyte count can be examined in peripheral blood film examination which is essential for hematology patients.

**2. General Component Required**

- Blood tubes with EDTA anticoagulant
- Clean glass slides
- Spreader with smooth edge
- Grease pencil
- Capillary tube
- Leishman's stain solution
- Distilled water or Sorensen's phosphate buffer (pH 6.8)




**Figure 1 :** Components of Blood Film Preparation

*Note: It is necessary to follow a standard procedure for blood collection to get the most accurate and trustworthy results of peripheral blood film examination.*

**3. Procedure**

**a. Peripheral Blood Film Preparation**

- 1) Gently invert the blood bottle for 5-7 times to mix thoroughly with the EDTA anticoagulant.
- 2) A small drop of blood is placed at one third of the slide with capillary tube.
- 3) The spreader or another clean slide is taken and placed it at an angle of 45 to the first slide.
- 4) The spreader is moved back until it meets the blood drop. When the blood has run along the edge of the spreader, the spreader is moved forward in a single rapid smooth movement pulling the blood after it.
- 5) The smear is allowed to dry in air and the blood film is ready for staining.



**Figure 2 :** Procedure of Blood Film Preparation



**Note:**

- The faster the spreader slide is moved, the longer and thinner the film will be. The slower the slide is moved, the shorter and thicker the slide will be
- An angle greater than 30° makes the smear thicker; less than 30° the smear is thinner.
- A small drop of blood may be insufficient to prepare a slide of sufficient length, too large a drop may cause the smear to extend beyond the length of the slide.

### **b. Blood Film Staining Method**


- 1) Fix the blood film by applying 3minutes 70% alcohol solution
- 2) Dry the film carefully
- 3) Dye the film using:

➤ **Quick Giemsa stain**

- 4) Place the prepared blood film in undiluted Giemsa Stain for 1 minutes.
- 5) Rinse in deionized water for 2 minutes.
- 6) Allow the slide to dry on absorbent paper and clean the back of the slide.

➤ **Leishman's stain**

- 4) The prepared thin blood film is fixed in Leishman's stain for 2minutes.
- 5) Twice volume of Sorensen's phosphate buffer or distilled water is added to the stain.
- 6) The film is washed and differentiated with Sorensen's phosphate buffer or D/W until the blood show salmon pink color.
- 7) It is drained and dried in air at RT. The back of the slide is cleaned and examined microscopically.




**Figure 3 :** Leishman Staining Method and Result Interpretation of Peripheral Blood Film

**Note:**

- Stain artifact such as debris and precipitates may be caused by over-staining (excess stain contact time) and inadequate washing under running water.
- Weak staining can be occurred because of using fatigued (old) stains. Stains should be checked before using and make new stain if required.

### **c. Blood film reading and result interpretation**

- 1) Clean the back of the and place the smear under the microscope.
- 2) Put a drop of oil.
- 3) Set the microscope the read the smear with the X100 oil immersion objective.
- 4) Observe the blood elements per
  - Type;
  - Size;
  - Shape;
  - Coloration;
  - Inclusion.



**Figure 4 :** Basic figure on result Interpretation on Peripheral Blood Film

**Note:**

- For the white blood cell differentiation, 100 white blood cells element should be observed
- A percentage of each white blood cell identified are reported in %
- The absolute value is calculated with: **WBC element% x WBC total count**



#### **4. Warnings and Precautions**

- Pre-analytical variables that can affect the quality of film must be controlled. Samples are best analyzed within 2 hours of blood collection. Delay in preparation of blood smear may allow for the degeneration of the cellular elements of blood and may lead to false result interpretation.
- Do not use the expired EDTA blood bottles.
- Do not perform the test if the sample contain blood clots.
- Do not use the unclean slide and spreader.
- Wear personal protective equipment, such as gloves and lab coats when handling patient sample. Handle all specimens as if they contain infectious agents.
- Do not smoke, drink, or eat while handling specimen. Long hair is to be either tied back or up. Placing pens and pencils in the mouth are avoided.
- Any spills of any biological sample of reagent to be decontaminated and disposed to minimize the risk of spreading infectious agents to the other person.
- Before leaving the laboratory, make ensure that your work area is clean and tidy as instructed.
- Wash hands thoroughly after the tests are done.
- Dispose of all specimens and materials used to perform the test as bio-hazard waste.
- Laboratory chemical and biohazard wastes must be handled and discarded in accordance with all local, state, and national regulations.

#### **References**

- University of Medical Technology University: Hematology General Practical Guidance Handbook
- Dace and Lewis: Practical Hematology, Twelfth Edition
- [http://www.blackwellpublishing.com/content/BPL/Images/Content\\_store/Sample\\_Chapter/1405142650/1405142650\\_4\\_001.pdf](http://www.blackwellpublishing.com/content/BPL/Images/Content_store/Sample_Chapter/1405142650/1405142650_4_001.pdf)