



GENERAL LABORATORY COMPONENTS REQUIRED IN HISTOPATHOLOGY

1. Introduction

Histopathology is the diagnosis and study of diseases of the tissues and involves examining tissues and/or cells under a microscope. Histopathologists are responsible for making tissue diagnoses and helping clinicians manage a patient’s care.

2. Items Used in Histology

 <p>Plastic (or) Camel hairbrush</p>	 <p>Coplin jar with lid</p>	 <p>Grease pencil</p>	 <p>Diamond pencil</p>
 <p>Timer</p>	 <p>MICROTOME KNIFE HONING GUIDE BLADE HANDLE Toe Heel</p> <p>Microtome knife with honing back (with or without handle)</p>	 <p>Plastic embedding ring with paraffin wax</p>	 <p>Stainless steel base mold (various sizes)</p>
 <p>Plastic tissue cassette</p>	 <p>Tissue block holder</p>	 <p>Single (or) Double slide mailer</p>	 <p>Microtome knife disposable blade</p>
 <p>Staining trough with lid</p>	 <p>Reagent bottle with glass lid or stopper</p>	 <p>Microtome</p>	 <p>Water bath</p>
 <p>H & E staining Reagent</p>	 <p>Clean glass slide and Cover slip</p>	 <p>DPX mounting media</p>	 <p>Microscope</p>

Table 1 Figures showing the items used in Histology



3. Common Fixatives used in Histology

The aim of fixation is to preserve cells or tissues in as near a life-like condition as possible, prevent autolysis and putrefaction, and protect the tissue from damage during subsequent processing. The most common used of fixatives are

- 1) 10% Formal saline
- 2) 10% Neutral buffered formalin

a. Preparation of **10% Formal saline**

Recommended for: Period of fixation: 24 hours or > depending on the tissue size. Indefinite period of fixation for large specimens with changing solutions every 3 months.

Formaldehyde, 40%	100 ml
Sodium chloride	8.5 g
Distilled water	900 ml
Formalin	10 parts
Formalin	90 arts

Procedure

- 8.5 g of NaCl is added to 500 ml DW in a beaker and stirred with a glass rod until it is completely dissolved. 100ml of formalin is added to NaCl solution.
- Then this mixture is transferred to 1L volumetric flask and made up to 1000 ml with DW.

b. Preparation of **10% Neutral buffered formalin (10% Buffered Formalin)**

Recommended for the preservation & storage of surgical, post-mortem and research specimens; the period of fixation is 24 hrs or longer.

Formaldehyde, 40%	100 ml
Sodium dihydrogen phosphate (anhydrous)	3.5g
Disodium hydrogen phosphate (anhydrous)	6.5g
Distilled water	900ml

Procedure

- Commercial formalin is buffered with 3.5g of acid sodium phosphate monohydrate and 6.5g of anhydrous disodium phosphate.
- This buffer formalin is transferred to the 1L volumetric flask and made up to 100ml with DW.

Reference

- <https://www.rcpath.org/discover-pathology/news/fact-sheets/histopathology.html>
- University of Medical Technology, Department of Medical Technology, Histopathology Practical Manual Handbook
- <https://www.radicalindia.com/histopathology-equipments.php>
- [https://en.wikipedia.org/wiki/Fixation_\(histology\)](https://en.wikipedia.org/wiki/Fixation_(histology))